RODAC PLATE

DESCRIZIONE

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BIBLIOGRAFIA

1. European Commission, Good Manufacturing Practices. Medical products for human and veterinary use. Ed. 1997.
2. APHA (1970). *Diagnostic Procedures and Reagents*.
3. Ronald M. Atlas. Handbook of Microbiological Media. 2nd Edition (1997).
4. Working Party on Control of Medicines and Inspections. Revision of the *Guide to Good Manufacturing Practice*. Annex on the *Manufacture of sterile medicinal products*. Brussels, 19/6/1995.
5. Project for ISO 18593. *Enumeration of aerobic bacteria using contact plates or dipslide and swabs methods*. Doc. ISO/TC34/SC 9 N 374, March, 1999.
6. EudraLex - Volume 4 Good manufacturing practice (GMP) Guidelines – Annex 1 (2008).

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RODAC PLATE

DESCRIPTION

RODAC plates are used for monitoring environmental hygiene (air and surfaces), to validate cleaning and disinfection operations, to apply Good Manufacturing Practices "GMP", to increase HACCP system and to train staff to hygiene. RODAC plates allow recovering and growth of microorganisms, which are present in the environment, on suitable culture medium.

CONTENT OF THE PACKAGES

- 1 instructions sheet
- 20 RODAC plates ready to use, packed one by one in blister packs of two pieces, five blister are wrapped in film thermally welded, and two packs of 10 plates contained in a cardboard box.

The irradiated packages (ref. 15340S, 15380S, 15336S, 15323S, 15365S, 15339S) show three envelopes: five blister are wrapped in film thermally welded, two packs are wrapped in vacuum bag. A cardboard box constitutes the final container of these envelopes. Success of irradiation is indicated by turning from yellow-orange to red of the colour of the circular indicator which is present on the label and on the two packs of five blister.

CONFIGURATIONS

RODAC plates proposed by Liofilchem are available in the following configurations:

RODAC plates without neutralizing

LEGIONELLA AGAR For <i>Legionella</i> spp. isolation	Ref. 15334	CETRIMIDE AGAR (acc. to Harmonized EP) For <i>Pseudomonas aeruginosa</i> isolation	Ref. 15332
LEGIONELLA AGAR (GVPC) For <i>Legionella</i> spp. isolation	Ref. 15376	ROSE BENGAL AGAR + CAF For yeasts and moulds isolation and enumeration	Ref. 15374
MAC CONKEY AGAR (acc. to Harmonized EP) For Gram-negatives isolation	Ref. 15329	R2A AGAR For heterotrophic microorganisms enumeration from drinking water	Ref. 15354
MANNITOL SALT AGAR (acc. to Harmonized EP) For staphylococci isolation	Ref. 15328	SABOURAUD DEXTROSE AGAR (acc. to Harmonized EP) For yeasts and moulds isolation	Ref. 15327
m-FAECAL COLIFORM AGAR For faecal coliforms isolation and enumeration	Ref. 15368	VIOLET RED BILE GLUCOSE AGAR (acc. to Harmonized EP / ISO 21528-2) For enterobacters isolation and enumeration	Ref. 15375
PLATE COUNT AGAR For total microbial count	Ref. 15325	VIOLET RED BILE LACTOSE AGAR For coliforms isolation and enumeration	Ref. 15326

Irradiated RODAC plates without neutralizing

TRYPTIC SOY AGAR CASEIN SOYA BEAN DIGEST AGAR (acc. to Harmonized EP) (Irradiated) For total microbial count	Ref. 15340S	SABOURAUD CAF AGAR (Irradiated) For yeasts and moulds isolation	Ref. 15380S
--	-------------	---	-------------

RODAC plates with neutralizing

BAIRD PARKER AGAR + NEUTRALIZING For <i>Staphylococcus aureus</i> isolation	Ref. 15331	SLANETZ B.A. + NEUTRALIZING * For faecal streptococci isolation and enumeration	Ref. 15362
PLATE COUNT AGAR + NEUTRALIZING For total microbial count on surfaces	Ref. 15336	T.S.A. + 0.1% PENASE + NEUTRALIZING For total microbial count	Ref. 15330
PLATE COUNT AGAR + TTC + NEUTRALIZING For total microbial count	Ref. 15360	V.R.B.G. + NEUTRALIZING For enterobacters isolation and enumeration	Ref. 15382
SABOURAUD CAF AGAR + NEUTRALIZING For yeasts and moulds isolation	Ref. 15365	V.R.B.L. + NEUTRALIZING For coliforms isolation and enumeration	Ref. 15361

Irradiated RODAC plates with neutralizing

PLATE COUNT AGAR + NEUTRALIZING (Irradiated) For total microbial count	Ref. 15336S	SABOURAUD CAF AGAR + NEUTRALIZING (Irradiated) For yeasts and moulds isolation	Ref. 15365S
SABOURAUD AGAR + NEUTRALIZING (Irradiated) * For yeasts and moulds isolation	Ref. 15323S	TRYPTIC SOY AGAR + NEUTRALIZING CASEIN SOYA BEAN DIGEST AGAR (acc. to Harmonized EP) (Irradiated) * For total microbial count	Ref. 15339S

PRINCIPLE OF THE METHOD

RODAC plates are used for monitoring environmental hygiene (air and surfaces), allowing recovering and growth of microorganisms, which are present in the environment, on suitable culture medium; beside they have the agar medium projecting over the opened plate, allowing a rapid contact of the medium with the surface under test. Disinfections and antiseptics, used for environmental sanitization, can inhibit microbial growth on the culture medium of the RODAC plates, so it is necessary that media contain substances inactivating disinfections and antiseptics: Lecithin, Histidine, Tween 80 and sodium thiosulphate neutralize disinfections and antiseptics commonly used such as ammonium quaternary compounds, phenolic derivatives, aldehydes, mercury compounds and active chlorine.

N.B.: For the products present in the previous tables and marked by the asterisk * (ref. 15323S, 15362, 15339S) neutralizing is constituted only by Lecithin and Tween 80.

TRYPTIC SOY AGAR CASEIN SOYA BEAN DIGEST AGAR (Irradiated), SABOURAUD CAF AGAR (Irradiated), PLATE COUNT AGAR + NEUTRALIZING (Irradiated), SABOURAUD DEXTROSE AGAR + NEUTRALIZING (Irradiated), SABOURAUD CAF AGAR + NEUTRALIZING (Irradiated) and TRYPTIC SOY AGAR+ NEUTRALIZING CASEIN SOYA BEAN DIGEST AGAR (Irradiated) plates are irradiated after packaging to assure the sterility of the outside envelopes (irradiation of the outside package is not performed for plates containing stainings).

TEST PROCEDURE

- For irradiated plates remove the first envelope from the package in the zone outside of the protected area. Remove the inside envelope in the protected area.
- Peel the blister and take the RODAC plate.
- Open the cap of the plate and gently push the agar surface on a level surface to examine for at least 10 seconds. Otherwise the plate can be used for air monitoring, using a "SAS" sampling and drawing a volume from 100 to 500 L.
- Cover the plate with the cap and incubate according to the conditions indicated in the following table:

	Aerobic mesophilic bacteria	Yeasts and moulds	Environmental bacteria	<i>Legionella</i> spp.
Time	72 hours	7 days	72 hours away from light and then for further 3 days at room temperature	Up to 10 days
Temperature	25-30 ± 1°C	25 ± 1°C	30 ± 1°C	36 ± 2°C
Atmosphere	aerobiosis	aerobiosis	aerobiosis	aerobiosis + humidity

INTERPRETATION OF RESULTS

Count Colony Forming Units developed on the agar surface manually or using automatic systems.

Evaluate the hygienic conditions of the microbiological control of surfaces using the following table:

Type of surface	Hygienic judgment	CFU/Plate
Floors, Walls, Work surfaces, Wash basins, Containers	Insufficient Acceptable Good	>50 26-50 0-25

Evaluate microbiological contamination of classified environment according to Annex 1 EU-GMP.

Recommended limits for microbial contamination		
Grade	Air sample CFU/m ³	Contact plates (surfaces) CFU/plate
A	<1	<1
B	10	5
C	100	25
D	200	50

QUALITY CONTROL

RODAC plates quality control provides for:

- Control of general characteristics, label and print
- Sterility control by incubation for 12 days at 25 ± 1°C and 12 days at 37 ± 1°C
- Microbiological control by inoculation of the plates with suspensions in physiological solution of the following strains:
Candida albicans ATCC 10231, *Enterococcus faecalis* ATCC 19433, *Escherichia coli* ATCC 25922, *Klebsiella pneumoniae* ATCC 13883, *Legionella pneumophila* ATCC 33152, *Proteus mirabilis* ATCC 25933, *Pseudomonas aeruginosa* ATCC 27853, *Saccharomyces cerevisiae* ATCC 9763, *Staphylococcus aureus* ATCC 25923. Inoculum for productivity 10²-10³ CFU/ml.

PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product must be used only by properly trained operators.

STORAGE

10-25°C (2-8°C for BAIRD PARKER AGAR + NEUTRALIZING ref. 15331 and LEGIONELLA AGAR (GVPC) ref. 15376) away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

REFERENCES

- European Commission, Good Manufacturing Practices. Medical products for human and veterinary use. Ed. 1997.
- APHA (1970). *Diagnostic Procedures and Reagents*.
- Ronald M. Atlas. Handbook of Microbiological Media. 2nd Edition (1997).
- Working Party on Control of Medicines and Inspections. Revision of the *Guide to Good Manufacturing Practice*. Annex on the *Manufacture of sterile medicinal products*. Brussels, 19/6/1995.
- Project for ISO 18593. *Enumeration of aerobic bacteria using contact plates or dipslide and swabs methods*. Doc. ISO/TC34/SC 9 N 374, March, 1999.
- Current EudraLex - Volume 4 Good manufacturing practice (GMP) Guidelines – Annex 1 (2008).

TABLE OF SYMBOLS

LOT Batch code	 Temperature limitation	 Manufacturer	 Contains sufficient for <n> tests
REF Catalogue number	 Keep away from heat	 Use by	 Caution, consult accompanying documents

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LIOFILCHEM® S.r.l.