Blood Agar Base

Non-selective medium for the isolation and cultivation of many pathogenic and non-pathogenic microorganisms.

DESCRIPTION

Blood Agar Base is a general purpose medium which may be enriched with blood or serum to support the growth of a wide variety of fastidious organisms of clinical significance.

With added blood, this medium is specially designed for the determination of typical hemolytic reactions which are important differentiating criteria for streptococci, staphylococci, and other organisms.

Blood Agar Base may be employed without adding blood as a nutrient agar or as a medium for the short-term maintenance of stock cultures.

TYPICAL FORMULA

	(g/l)
Tryptose	10.0
Meat Extract	10.0
Sodium Chloride	5.0
Agar	15.0
Final pH 7.3 ± 0.2 at 25°C	

METHOD PRINCIPLE

Tryptose and meat extract provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Sodium chloride maintains the osmotic balance of the medium. Agar is the solidifying agent.

PREPARATION

<u>Dehydrated medium</u>	Suspend 40 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes.
<u>Medium in bottles</u>	Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes.

For blood agar, add 5-10% sterile defibrinated sheep (ref. 83394) or horse (ref. 83294) blood to the Base cooled at 45-50°C after sterilization. Mix with gentle rotation and pour into petri dishes or other containers.

Boiled blood agar (chocolate agar) is prepared by heating after adding the blood. Heat the culture medium for about 10 minutes at approx. 80°C with frequent swirling until it turns brownish. Chocolate agar supplemented with 10% defibrinated blood is suited for isolating *Haemophilus* and *Neisseria* species.

For the selective isolation of tubercle bacilli the addition of 1% glycerol and 25% human blood is recommended.

TEST PROCEDURE

Process each specimen as appropriate, and inoculate over the agar surface. Streak for isolation with an inoculating loop, then stab the agar several times to deposit beta-hemolytic streptococci beneath the agar surface. Subsurface growth will display the most reliable hemolytic reactions demonstrating both oxygen-stable and oxygen-labile streptolysins. Incubate plates aerobically, anaerobically or under CO₂ enriched atmosphere in accordance with established laboratory procedures.

INTERPRETING RESULTS

Examine plates for growth and hemolytic reactions after 18-48 hours of incubation. Four different types of hemolysis on blood agar media can be described:

1. Alpha (α)-hemolysis is the reduction of hemoglobin to methemoglobin in the medium surrounding the colony, causing a greenish discolorization of the medium.
2. Beta (β)-hemolysis is the lysis of red blood cells, resulting in a clear zone surrounding the colony.
3. Gamma (γ)-hemolysis indicates no hemolysis. No destruction of red blood cells occurs, and there is no change in the medium.
4. Alpha-prime (α')-hemolysis is a small zone of complete hemolysis that is surrounded by an area of partial lysis.

The haemolytic reactions of organisms inoculated on to this medium will be affected by the animal blood used e.g. horse or sheep and the incubation conditions e.g. aerobic, capnoeic or anaerobic. When horse blood is added to the medium *Haemophilus haemolyticus* colonies will produce beta-haemolysis and mimic *Streptococcus pyogenes*.

APPEARANCE

Dehydrated medium: free-flowing, homogeneous, light beige.

Prepared medium: before adding blood, clear and yellowish-brown, then cherry red.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store bottles at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.

Medium in bottles: 2 years.

QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 10-100 CFU.

Incubation conditions: 35 ± 2°C for 18-24 hours in aerobic atmosphere supplemented with CO₂.

QC Table.

Microorganism		Growth	Hemolysis
<i>Streptococcus pyogenes</i>	ATCC® 19615	Good	Beta
<i>Streptococcus pneumoniae</i>	ATCC® 6305	Good	Alpha
<i>Staphylococcus aureus</i>	ATCC® 25923	Good	Beta
<i>Escherichia coli</i>	ATCC® 25922	Good	---

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. Quality Control for Commercially Prepared Microbiological Media (2004) - 3rd ed. M22-A3. Clinical and Laboratory Standards Institute - CLSI (NCCLS), Wayne, PA.
2. Ruoff, Whiley and Beighton. (1999). In Murray, Baron, Pfaller, Tenover and Tenover (ed.), Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.
3. Isenberg (ed).1992. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
4. Facklam R. R. (1980) in Manual of Clinical Microbiology. Eds. Lennette E. H., Balows A., Hausler W. J. & Truant J. P. 3rd Edn. Amer. Soc. for Microbiology. Washington D.C. pp. 88-110.
5. M.S. Tarshis, A.W. Frisch (1951) Blood media for the cultivation of Mycobacterium tuberculosis, Amer. J. Clin. Pathol., 21:101.

PRESENTATION

PRESENTATION		Contents	Ref.
Blood Agar Base	Bottles	6 x 500 ml bottles	470130
Blood Agar Base	Bottles	6 x 100 ml bottles	402140
Blood Agar Base	Dehydrated medium	500 g of powder	610005
Blood Agar Base	Dehydrated medium	100 g of powder	620005

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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BIBLIOGRAFIA

1. Quality Control for Commercially Prepared Microbiological Media (2004) - 3rd ed. M22-A3. Clinical and Laboratory Standards Institute - CLSI (NCCLS), Wayne, PA.
2. Ruoff, Whiley and Beighton. (1999). In Murray, Baron, Pfaller, Tenover and Tenover (ed.), Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.
3. Isenberg (ed).1992. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
4. Facklam R. R. (1980) in Manual of Clinical Microbiology. Eds. Lennette E. H., Balows A., Hausler W. J. & Tenover J. P. 3rd Edn. Amer. Soc. for Microbiology. Washington D.C. pp. 88-110.
5. M.S. Tarshis, A.W. Frisch (1951) Blood media for the cultivation of Mycobacterium tuberculosis, Amer. J. Clin. Pathol., 21:101.

PRESENTAZIONE		Contenuto	Ref.
Blood Agar Base	Flaconi	Flaconi 6 x 500 ml	470130
Blood Agar Base	Flaconi	Flaconi 6 x 100 ml	402140
Blood Agar Base	Terreno disidratato	500 g di polvere	610005
Blood Agar Base	Terreno disidratato	100 g di polvere	620005

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- Quality Control for Commercially Prepared Microbiological Media (2004) - 3rd ed. M22-A3. Clinical and Laboratory Standards Institute - CLSI (NCCLS), Wayne, PA.
- Ruoff, Whiley and Beighton. (1999). In Murray, Baron, Pfaller, Tenover and Tenover (ed.), Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.
- Isenberg (ed). 1992. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
- Facklam R. R. (1980) in Manual of Clinical Microbiology. Eds. Lennette E. H., Balows A., Hausler W. J. & Tenover J. P. 3rd Edn. Amer. Soc. for Microbiology. Washington D.C. pp. 88-110.
- M.S. Tarshis, A.W. Frisch (1951) Blood media for the cultivation of Mycobacterium tuberculosis, Amer. J. Clin. Pathol., 21:101.

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		Contenido	Ref.
Blood Agar Base	Botellas	6 x 500 ml botellas	470130
Blood Agar Base	Botellas	6 x 100 ml botellas	402140
Blood Agar Base	Medio deshidratado	500 g de polvo deshidratado	610005
Blood Agar Base	Medio deshidratado	100 g de polvo deshidratado	620005

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