Mueller Hinton Broth

Medium for antimicrobial susceptibility testing of aerobic organisms.

TYPICAL FORMULA	(g/l)
Beef Extract	2.0
Acid Hydrolysate of Casein	17.5
Starch	1.5
Final pH 7.3 ± 0.1 at 25°C	

DESCRIPTION

Mueller Hinton Broth is used for the cultivation of a wide variety of fastidious and nonfastidious microorganisms. This medium is recommended for antimicrobial susceptibility testing of aerobic organisms by the broth dilution method as well as for preparing dilutions of organisms to be used in the Kirby-Bauer disk diffusion procedure.

PRINCIPLE

Casamino acids and beef extract provide nutrients required for the organisms growth. Agar is the solidifying agent. Starch acts as a colloid that protects against toxic material in the media

PREPARATION

Suspend 21.0 g of powder in 1 liter of distilled or deionized water. Heat to dissolve completely. Dispense into final containers. Sterilize in autoclave at 121 °C for 15 minutes.

TECHNIQUE

Organisms to be subcultured must first be isolated in pure culture on an appropriate solid medium. Transfer growth from the isolation medium to Mueller Hinton Broth using standard bacteriologic techniques.

For enrichment purposes, inoculate the specimen onto primary media and then into the broth according to recommended procedures. Incubate tubes and bottles at 35°C under conditions appropriate for the organism being cultured.

INTERPRETATION OF RESULTS

Growth in broth media is indicated by the presence of turbidity compared with an uninoculated control.

Note: Many factors can affect results: inoculum size, rate of growth, medium formulation and pH, length of incubation and incubation environment, disk content and drug diffusion rate, and measurements of endpoints. Therefore, strict adherence to protocol is required to ensure reliable results. A pH outside the expected range may adversely affect the results of the susceptibility test. If the pH is too low, aminoglycosides and macrolides will appear to lose potency, others may appear to have excessive activity. The opposite effects are possible if the pH is too high.

STORAGE AND TRANSPORT CONDITIONS

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident. Store prepared plates at 2-8°C away from light.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for professional use only and must be used by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

- Wood, G.L. and J.A. Washington (1995) Antimicrobial Susceptibility Tests; dilution and disk diffusion methods. In Murray, P.R., et al. *Manual of Clinical Microbiology*, 7th ed. American Society for Microbiology, Washington, D.C.; 1327-1341.
- Tille, P., et al. *Bailey and Scott's Diagnostic Microbiology*, C.V. Mosby Company, St. Louis, MO.
- Barry and Fay (1973) *Am. J. Clin. Pathol.*; 50:196.
- Bauer, A.W., W.M.M. Kirby, et al. (1966) *Am. J. Clin. Pathol.*; 45:493-496.
- Isenberg, H.D. *Clinical Microbiology Procedures Handbook*, Vol. I, II & III. American Society for Microbiology, Washington, D.C.
- Koneman, E.W., et al. *Color Atlas and Textbook of Diagnostic Microbiology*, J.B. Lippincott Company, Philadelphia, PA.
- *Methods for Dilution Antimicrobial Test for Bacteria that Grow Aerobically*, 4th ed., M7-A5. 1997. Clinical Laboratory Standards Institute (CLSI - formerly NCCLS), Villanova, PA.
- Mueller, J.H. and J. Hinton (1941) A protein-free medium for primary isolation of the *Gonococcus* and *Meningococcus*. *Proc. Soc. Exp. Biol. and Med.*; 48:330-333.
- Standard Disk Susceptibility Test, The Federal Register, September 30, 1972; 37(191):20527-20529.

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10-30°C

PACKAGING

Ref.	Content	Packaging
610034	500 g	500 g of powder in plastic bottle
620034	100 g	100 g of powder in plastic bottle

pH OF THE MEDIUM

7.3 ± 0.1

USE

Mueller Hinton Broth is used for the cultivation of a wide variety of microorganisms and for antimicrobial susceptibility testing by broth dilution and disk diffusion methods

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Powder medium

Appearance: free-flowing, homogeneous

Colour: beige

Ready-to-use medium

Appearance: slightly opalescent

Colour: light to medium amber

SHELF LIFE

4 years

QUALITY CONTROL

- Control of general characteristics, label and print
- Microbiological control
Inoculum for productivity: ≤100 CFU
Incubation Conditions: 18-24 h at 36 ± 1°C

Microorganism		Growth
<i>Escherichia coli</i>	ATCC® 25922	Good
<i>Enterococcus faecalis</i>	ATCC® 29212	Good
<i>Staphylococcus aureus</i>	ATCC® 25923	Good

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult instructions for use	 Do not reuse

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LIOFILCHEM® S.r.l.