MacConkey Sorbitol Agar

Selective and differential medium for detection of *E. coli* O157 from clinical samples and other materials.

DESCRIPTION

MacConkey Sorbitol Agar is a partially selective and differential medium used for the isolation of enterohemorrhagic *Escherichia coli* O157:H7, from clinical specimens, food and environmental samples.

TYPICAL FORMULA

	(g/l)
Pancreatic Digest of Gelatin	17.0
Peptone from Meat	1.5
Peptone from Casein	1.5
D-Sorbitol	10.0
Sodium Chloride	5.0
Bile Salts	1.5
Agar	15.0
Neutral Red	0.03
Crystal Violet	0.001
Final pH 7.1 ± 0.2 at 25°C	

METHOD PRINCIPLE

Pancreatic digest of gelatin and peptones from meat and casein provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Sorbitol is the fermentable carbohydrate. Sodium Chloride maintains the osmotic balance of the medium. Bile salts and crystal violet are the selective agents, inhibiting Gram-positive organisms and allowing Gram-negative bacteria to grow. Agar is the solidifying agent. Neutral red is the pH indicator.

PREPARATION

- Dehydrated medium Suspend 51.5 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil for 1 minute shaking frequently until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes.
- Medium in bottles Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes.

TEST PROCEDURE

Inoculate the plates by directly streaking the specimen on the agar surface or spread the sample from an enrichment culture. Incubate aerobically at 35 ± 2°C for 18-24 h

INTERPRETING RESULTS

Sorbitol-fermenting organisms, such as common *E. coli* strains, produce pink to red colonies.

Organisms that do not ferment sorbitol, such as *E. coli* O157:H7, grow as colorless colonies. Serological or molecular tests are required for confirmation of the serotype and/or of the toxins.

Enterococci, Staphylococci and other Gram-positive bacteria are partially or completely inhibited.

APPEARANCE OF THE MEDIUM

Dehydrated medium: free-flowing, homogeneous, beige-pink.

Prepared medium: slightly opalescent, pinkish-red.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store bottles and prepared plates at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.
 Medium in bottles: 2 years.
 Ready-to-use plates: 6 months.

QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC table.
 Inoculum for productivity: 50-100 CFU
 Inoculum for selectivity: 10⁴-10⁶ CFU
 Incubation conditions: aerobically at 35 ± 2°C for 18-24 hours.

QC Table.

Microorganism		Growth	Specification
<i>Escherichia coli</i> O157:H7 (sorbitol-negative)	ATCC® 33150	Good	Colorless colonies
<i>Escherichia coli</i> (sorbitol-positive)	ATCC® 25922	Good	Pink colonies
<i>Enterococcus faecalis</i>	ATCC® 29212	Partially to completely inhibited	Very small, opaque colonies

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* diagnostic use only and must be used by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

- Murray, Baron, Jorgensen, Landry and Pfaller ed. (2007) Manual of clinical microbiology, 9th ed. American Society for Microbiology, Washington, D.C.
- Kaper, J.B., O'Brien, A.D. ed. (1998) *Escherichia coli* O157:H7 and Other Shiga Toxin- Producing *E. coli* strains. American Society for Microbiology, Washington, DC, USA.
- Mortlock, S (1994) Recovery of *Escherichia coli* O157:H7 from mixed suspensions: evaluation and comparison of pre-coated immunomagnetic beads and direct plating. *Brit. J. Biomed. Sci.* 51: 207-214.
- MacConkey A. (1905) Lactose-fermenting bacteria in faeces. *J. Hygiene* 8:333-379.

PRESENTATION

		Contents	Ref.
MacConkey Sorbitol Agar	90 mm ready-to-use plates	20 plates	10005
MacConkey Sorbitol Agar	90 mm ready-to-use plates	100 plates	10005*
MacConkey Sorbitol Agar	Bottles	6 x 100 ml bottles	402450
MacConkey Sorbitol Agar	Dehydrated medium	500 g of powder	610108
MacConkey Sorbitol Agar	Dehydrated medium	100 g of powder	620108

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Medical Diagnostic Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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