Corn Meal Agar

Medium for the cultivation of fungal species and demonstration of chlamyospore production.

DESCRIPTION

Corn Meal Agar is a medium used for the cultivation of fungi and identification of *C. albicans* based on the formation of chlamyospores.

TYPICAL FORMULA

	(g/l)
Corn Meal Infusion (from 50 g whole maize)	2.0
Agar	15.0
Final pH 6.0 ± 0.2 at 25°C	

METHOD PRINCIPLE

Infusion from corn meal provides amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Agar is the solidifying agent.

PREPARATION

Dehydrated medium Suspend 17 g of the powder in 1 liter of distilled or deionized water (*). Mix well. Heat to boil shaking frequently until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes.

(*). The addition of 1% polysorbate 80 (Tween 80) may enhance chlamyospore formation by *Candida albicans*. With the addition of 1% dextrose a more luxuriant growth of some fungi can be seen.

TEST PROCEDURE

Using a sterile inoculating needle or loop, touch the yeast colony and immediately make a deep cut in the agar. Place a sterile coverslip over the line of inoculum. Incubate aerobically at 25 ± 2°C for 2-7 days.

Notice that for chlamyospore production some *Candida* spp may take up to 2 weeks.

INTERPRETING RESULTS

Using a low power objective examine daily through the coverslip for growth, pigmentation and chlamyospore formation.

APPEARANCE

Dehydrated medium: free-flowing, homogeneous, beige to yellow.
Prepared medium: trace to moderately hazy, grey to yellow-white.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store prepared plates at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.
Ready-to-use plates: 6 months.

QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 50-100 CFU

Incubation conditions: aerobically at 25 ± 2°C for 2-7 days.

QC Table.

Microorganism		Growth	Chlamydo spores
<i>Aspergillus brasiliensis</i>	ATCC® 16404	Good	N/A
<i>Candida albicans</i>	ATCC® 10231	Good	Positive
<i>Saccharomyces cerevisiae</i>	ATCC® 9763	Good	Negative

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

- Warren, N., and H. J. Shadomy (1995) Yeasts of medical importance, p. 617-629. In P.R. Murray, E. J. Baron, M.A. Pfaller, F. C. Tenover, and R. H. Tenkens (eds.). Manual of clinical microbiology, 6th ed., American Society for Microbiology, Washington, D.C.
- Baron, E. J. and S. M. Finegold (1990) Formulas and preparation of culture media and reagents. Bailey & Scott's Diagnostic Microbiology, 8th ed. The C. V. Mosby Company, St. Louis, MO.
- Duncan, J., and J. Floeder (1963) A Comparison of media for the production of chlamydo spores by *Candida albicans*. Am. J. Med. Tech. 29:199-206.
- Kelly, J. P., and F. Funigeillo (1959) *Candida albicans*: A study of media designed to promote chlamydo spore production. J. Lab & Clin. Med. 53:807-809.

PRESENTATION

		Contents	Ref.
Corn Meal Agar	90 mm ready-to-use plates	20 plates	11507
Corn Meal Agar	Dehydrated medium	500 g of powder	610123
Corn Meal Agar	Dehydrated medium	100 g of powder	620123

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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BIBLIOGRAFIA

- Warren, N., and H. J. Shadomy (1995) Yeasts of medical importance, p. 617-629. In P.R. Murray, E. J. Baron, M.A. Pfaller, F. C. Tenover, and R. H. Tenkels (eds.). Manual of clinical microbiology, 6th ed., American Society for Microbiology, Washington, D.C.
- Baron, E. J. and S. M. Finegold (1990) Formulas and preparation of culture media and reagents. Bailey & Scott's Diagnostic Microbiology, 8th ed. The C. V. Mosby Company, St. Louis, MO.
- Duncan, J., and J. Floeder (1963) A Comparison of media for the production of chlamydospores by *Candida albicans*. Am. J. Med. Tech. 29:199-206.
- Kelly, J. P., and F. Funigeillo (1959) *Candida albicans*: A study of media designed to promote chlamydospore production. J. Lab & Clin. Med. 53:807-809.

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