Campylobacter Karmali Agar Base

Basal medium for the isolation of *Campylobacter jejuni* and *Campylobacter coli* from clinical and nonclinical specimens.

TYPICAL FORMULA	(g/l)
Peptospecial	23.0
Starch	1.0
Activated Charcoal	4.0
Sodium Chloride	5.0
Agar	15.0
Final pH 7.4 ± 0.2 at 25°C	

DESCRIPTION

Campylobacter Karmali Agar Base is a blood-free, charcoal-based medium used with supplements for the selective isolation of *Campylobacter jejuni* and *Campylobacter coli* from feces and other materials.

PRINCIPLE

Peptospecial provides amino acids, nitrogen, carbon, minerals, vitamins and other nutrients which support the growth of microorganisms. Corn starch is a source of energy. Activated charcoal helps to neutralize the toxic metabolic products formed in the medium. Sodium chloride maintains the osmotic balance of the medium. Agar is the solidifying agent.

Supplementation with Campylobacter Karmali Supplement (ref. 81036), containing cefoperazone, vancomycin, cycloheximide, hemin and sodium pyruvate, inhibits the contaminant microbial flora in fecal specimens improving, at the same time, the recovery of *Campylobacter* spp.

PREPARATION

Suspend 48.0 g of powder in 1 liter of deionized or distilled water. Bring to boil and shake until completely dissolved. Sterilize at 121°C for 15 minutes. Cool up to 45-50°C. Aseptically, add rehydrated content of 2 vials (10 ml) of Campylobacter Karmali Supplement. Pour in Petri dishes.

TECHNIQUE

Inoculate the plates by streaking directly the sample onto the agar surface. Incubate under microaerobic atmosphere (containing approximately 5% O₂, 10% CO₂ and 85% N₂) at 42°C for a minimum of 48 hours.

INTERPRETATION OF RESULTS

The typical colonies of *Campylobacter* are small, greyish and flat with irregular edges. Colonies have a tendency to spread, especially when initially isolated from fresh clinical specimens.

STORAGE AND TRANSPORT CONDITIONS

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident. Store prepared plates at 2-8°C away from light.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *in vitro* diagnostic use only and must be used by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

- Vanderzant C., D.F. Splittstoesser (1992) Compendium of Methods for the Microbiological Examination of Foods, 3rd Ed. APHA, Washington D.C.
- Karmali M.A., A.E. Simor, M. Roscoe, P.C. Fleming, S.S. Smith, J. Lane (1986) J. Clin. Micro. 23:456-459.
- Skirrow M.B., J. Benjamin (1980) J. Clin. Path. 33:1122.

PRODUCT SPECIFICATIONS

NAME

Campylobacter Karmali Agar Base

PRESENTATION

Dehydrated medium

STORAGE

10-30°C

PACKAGING

Ref.	Content	Packaging
610200	500 g	500 g of powder in plastic bottle
620200	100 g	100 g of powder in plastic bottle

pH OF THE MEDIUM

7.4 ± 0.2

USE

Campylobacter Karmali Agar Base is a blood-free, charcoal-based medium used with supplements for the selective isolation of *Campylobacter jejuni* and *Campylobacter coli* from clinical and nonclinical specimens

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Powder medium

Appearance: free-flowing, homogeneous

Colour: beige

Ready-to-use medium

Appearance: opaque

Colour: black

SHELF LIFE

4 years

QUALITY CONTROL

- Control of general characteristics, label and print
- Microbiological control
Inoculum for productivity: 50-100 CFU
Inoculum for selectivity: 10⁴-10⁶ CFU
Incubation Conditions: 24-48 h at 42°C, under microaerobic atmosphere

Microorganism

Campylobacter jejuni

ATCC® 33291

Growth

Good

Escherichia coli

ATCC® 25922

Inhibited

TABLE OF SYMBOLS

 LOT	Batch code	 IVD	<i>In vitro</i> Diagnostic Medical Device		Manufacturer		Use by		Fragile, handle with care
 REF	Catalogue number		Temperature limitation		Contains sufficient for <n> tests		Caution, consult instructions for use		Do not reuse

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LIOFILCHEM® S.r.l.