Cooked Meat Medium

Medium for the cultivation of anaerobes,
from clinical and nonclinical samples.

DESCRIPTION

Cooked Meat Medium is used for the cultivation and maintenance of anaerobic organisms, especially pathogenic clostridia, differentiating saccharolytic from proteolytic species.

This medium permits growth from a small inoculum, important for clinical specimens. It is also recommended in standard methods for food testing for the enumeration and identification of *Clostridium perfringens*.

TYPICAL FORMULA

	(g/l)
Heart Tissue Granules	98.0
Peptic Digest of Animal Tissue	20.0
Dextrose	2.0
Sodium Chloride	5.0

Final pH 7.2 ± 0.2 at 25°C

METHOD PRINCIPLE

Heart tissue granules and peptic digest of animal tissue provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Solid meat particles provide favorable growth conditions for anaerobes due to reducing action of sulfhydryl (-SH) groups of muscle protein. Sulfhydryl groups are more accessible in denatured proteins, therefore the meat particles are cooked for use in the medium. Dextrose is the fermentable carbohydrate. Sodium chloride maintains the osmotic balance of the medium.

PREPARATION

<u>Dehydrated medium</u>	Place 1.25 g of the particles in a test tube and add 10 milliliters of distilled or deionized water. Sterilize in autoclave at 121°C for 15 minutes.
<u>Medium in tubes</u>	To obtain proper incubation conditions for growth of strict anaerobes, Cooked Meat Medium should be reduced prior to use by boiling the tubes for 10 min with caps loosened and cooling with tightened caps to room temperature before inoculation. Alternatively, an efficient and easy way to obtain suitable anaerobic conditions is through the use of the GasPak EZ anaerobic system (ref. 96726, 96724 or 96728).

TEST PROCEDURE

Inoculate the tube in the area of meat particles (bottom of the tube) with isolate from a fresh subculture medium. Incubate anaerobically at 35 ± 2°C for up to 7 days.

INTERPRETING RESULTS

Turbidity of the medium indicates microbial growth. In the cultivation of clostridia, saccharolytic organisms usually produce acid and gas. Growth of proteolytic organisms is generally characterized by blackening and dissolution of the meat particles.

APPEARANCE

Dehydrated medium: brown homogeneous pellets.

Prepared medium: amber, clear to very slightly opalescent supernatant fluid over insoluble pellets.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store tubes at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.

Medium in tubes: 2 years.

QUALITY CONTROL

Tubes are inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 10^2 - 10^3 CFU.

Incubation conditions: anaerobically at $35 \pm 2^\circ\text{C}$ for 40-48 hours.

QC Table.

Microorganism		Growth
<i>Clostridium perfringens</i>	ATCC® 12924	Good
<i>Clostridium sporogenes</i>	ATCC® 11437	Good
<i>Staphylococcus aureus</i>	ATCC® 25923	Good

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. U.S. Food and Drug Administration (2001) Bacteriological analytical manual, online. AOAC, International, Gaithersburg, Md.
2. Forbes B.A. and P.A. Granato (1995) Processing specimens for bacteria, p. 265-281. In P.R. Murray, M.A. Tenover and R.H. Tenover (ed.) Manual of clinical microbiology 6th ed. American Society for Microbiology, Washington, D.C.
3. Robertson M (1916) Notes upon certain anaerobes isolated from wounds. J. Pathol. Bacteriol. 20:327-349.

PRESENTATION

PRESENTATION		Contents	Ref.
Cooked Meat Medium	Tubes	20 x 20 ml tubes	24071
Cooked Meat Medium	Dehydrated medium	500 g of powder	610372

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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BIBLIOGRAFIA

1. I.U.S. Food and Drug Administration (2001) Bacteriological analytical manual, online. AOAC, International, Gaithersburg, Md.
2. Forbes B.A. and P.A. Granato (1995) Processing specimens for bacteria, p. 265-281. In P.R. Murray, M.A. Tenover and R.H. Tenover (ed.) Manual of clinical microbiology 6th ed. American Society for Microbiology, Washington, D.C.
3. Robertson M (1916) Notes upon certain anaerobes isolated from wounds. J. Pathol. Bacteriol. 20:327-349.

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