## Chromatic™ Staph aureus

Chromogenic medium for the isolation and differentiation of *Staphylococcus aureus* directly from clinical and nonclinical specimens.

TYPICAL FORMULA	(g/l)
Peptone and Yeast Extract	40.0
Sodium Chloride	10.0
Chromogenic Mix	7.0
Agar	15.0
Final pH 7.0 ± 0.2 at 25°C	

### DESCRIPTION

Chromatic™ Staph aureus is a selective chromogenic medium used with supplements for the isolation and differentiation of *S. aureus* directly from clinical and nonclinical specimens.

### PRINCIPLE

Peptone and yeast extract supply amino acids, nitrogen, carbon, minerals, vitamins and other nutrients which support the growth of microorganism. Sodium chloride maintains the osmotic balance of the medium enhancing growth of *S. aureus* as well. Chromogenic mix allows the identification of microorganisms on the basis of the colony color and morphology. Agar is the solidifying agent.

Supplementation with Chromatic™ Staph aureus Supplement (ref. 81085) inhibits the growth of most organisms other than staphylococci enhancing bacterial differentiation as well.

### PREPARATION

Suspend 72.0 g of powder in 1 liter of deionized or distilled water. Bring to boil and shake until completely dissolved. Sterilize at 121°C for 15 minutes. Cool up to 45-50°C. Aseptically, add rehydrated content of 2 vials (10 ml) of Chromatic™ Staph aureus Supplement (ref. 81085). Pour in Petri dishes.

### TECHNIQUE

Inoculate the plates by streaking directly the sample onto the agar surface. Incubate aerobically at 37°C for 18-24 hours.

### INTERPRETATION OF RESULTS

*S. aureus* produces pink-purplish colonies while other staphylococci cultivate with greenish colonies.

### STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident. Store prepared plates at 2-8°C away from light.

### WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *in vitro* diagnostic use only and must be used by properly trained operators.

### DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

### REFERENCES

1. O. Gaillot. Evaluation of CHROMagar Staph aureus, a new chromogenic medium, for isolation and presumptive identification of *Staphylococcus aureus* from human clinical specimens. Laboratoire de Bacteriologie-Virologie, Hopital Necker-Paris, France, rapport de Septembre 1999.
2. Carricao A., A.M. Freydiere et al. Performance of the Chromogenic medium CHROMagar Staph aureus and the Staphylochrom Coagulase Test in the detection and identification of *Staphylococcus aureus* in clinical specimens. Journal of clinical microbiology, July 2001, P.2581-2583, Vol:39 N°7.

## PRODUCT SPECIFICATIONS

### NAME

Chromatic™ Staph aureus

### PRESENTATION

Dehydrated medium

### STORAGE

10-30°C

### PACKAGING

Ref.	Content	Packaging
610616	500 g	500 g of powder in plastic bottle
620616	100 g	100 g of powder in plastic bottle

### pH OF THE MEDIUM

7.0 ± 0.2

### USE

Chromatic™ Staph aureus is a chromogenic medium used with supplements for the isolation and differentiation of *S. aureus* directly from clinical and nonclinical specimens

### TECHNIQUE

Refer to technical sheet of the product

### APPEARANCE OF THE MEDIUM

#### Powder medium

Appearance: free-flowing, homogeneous

Colour: beige

#### Ready-to-use medium

Appearance: slightly opalescent

Colour: light amber

### SHELF LIFE

2 years

### QUALITY CONTROL

- Control of general characteristics, label and print
- Microbiological control  
Inoculum for productivity: 10-100 CFU/ml  
Inoculum for selectivity: 10<sup>4</sup>-10<sup>5</sup> CFU/ml  
Inoculum for specificity: ≤10<sup>4</sup> CFU/ml  
Incubation Conditions: 18-24 h at 35 ± 2°C, in aerobiosis

Microorganism		Growth	Colony color
<i>Staphylococcus aureus</i>	ATCC® 6538	Good	Pink-purplish
<i>Staphylococcus sciuri</i>	ATCC® 29062	Good	Greenish
<i>Staphylococcus xylosus</i>	ATCC® 29971	Good	Greenish
<i>Escherichia coli</i>	ATCC® 25922	Inhibited	---

### TABLE OF SYMBOLS

 <b>LOT</b>	Batch code	 <b>IVD</b>	<i>In vitro</i> Diagnostic Medical Device		Manufacturer		Use by		Fragile, handle with care
 <b>REF</b>	Catalogue number		Temperature limitation		Contains sufficient for <n> tests		Caution, consult instructions for use		Do not reuse

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**LIOFILCHEM® S.r.l.**