## Chromatic™ *Bacillus cereus*

Chromogenic medium for detection of *Bacillus cereus* in food and environmental samples.

TYPICAL FORMULA	(g/l)
Peptone	10.0
Yeast Extract	4.0
Di-Sodium Hydrogen Phosphate	2.5
Potassium Dihydrogen Phosphate	0.3
Sodium Pyruvate	4.0
Chromogenic Mix	3.2
Agar	15.0
Final pH 7.2 ± 0.2 at 25°C	

### DESCRIPTION

Chromatic™ *Bacillus cereus* is a chromogenic medium used with supplements for the isolation, differentiation and enumeration of *Bacillus cereus* from a wide variety of foods including both raw and cooked food as well as sauces, soups and other food mixtures.

### PRINCIPLE

Peptone is a source of amino acids, nitrogen, carbon, minerals and vitamins essential for bacterial growth. Yeast extract provides vitamins particularly of B-group. Di-sodium phosphate and potassium phosphate are the buffer system. Sodium pyruvate is an energy source. The chromogenic mix allows the identification of bacteria on the basis of the colony color. Agar is the solidifying agent.

This medium is completed after the addition of Chromatic™ *Bacillus cereus* Supplement (ref. 81097) consisting of an enrichment supplement and a selective supplement. Selectivity is due to the combination of two antibiotics: polymyxin B and trimethoprim. Polymyxin B inhibits most Gram-negative organisms and some Gram-positive bacteria including some bacilli other than *B. cereus*. Trimethoprim, which blocks folic acid synthesis necessary for cellular replication, is active against many Gram-positive bacteria including *S. aureus*, *Enterococcus* spp and some non-*Bacillus cereus* spp.

### PREPARATION

Suspend 39 g of powder in 1 liter of deionized or distilled water. Bring to boil and shake until completely dissolved. Sterilize at 121°C for 15 minutes. Cool up to 45-50°C. Aseptically, add rehydrated content of 2 vials (10 ml) of Selective Supplement and 2 vials (40 ml) of Enrichment Supplement, both supplements contained in Chromatic™ *Bacillus cereus* Supplement. Mix well and pour in Petri dishes.

### TECHNIQUE

Prepare the sample by diluting 1:10 (w/v) with Buffered Peptone Water (ref. 24099) or Ring's Solution (ref. 81059). Homogenize for 1 minute by using an appropriate laboratory blender and spread the material over the agar surface. Alternatively, inoculate the medium by direct streaking of the sample onto plate. Incubate aerobically at 30°C for 18-24 hours.

### INTERPRETATION OF RESULTS

Examine for typical blue-green colonies of *Bacillus cereus* surrounded by a opaque halo. Confirm the presumptive identification by a validated method.

Note that *Bacillus cereus* strains which do not produce the enzyme  $\beta$ -glucosidase cannot cleave the chromogenic substrate thereby forming white colonies with opaque halo.

### STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident. Store prepared plates at 2-8°C away from light.

### WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for professional use only and must be used by properly trained operators.

### DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

### REFERENCES

1. ISO 7932:2004. Horizontal method for the enumeration of presumptive *Bacillus cereus*.
2. ISO 16140:2003. Microbiology of food and animal feeding stuffs - Protocol for the validation of alternative
3. Foodborne Pathogenic Microorganisms and Natural Toxins Handbook *Bacillus cereus* and other *Bacillus* spp. (2003) U.S. Food & Drug Administration (C.F.S.A.N)
4. Handbook of Culture Media for Food Microbiology (2003) Volume 37. Chapter 4. Media for *Bacillus* spp. and related genera relevant to foods. Edited by Corry, J. E. L., Curtis, G. D. W. and Baird, R. M. Publisher - Elsevier, Amsterdam.

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10-30°C

### PACKAGING

Ref.	Content	Packaging
610628	500 g	500 g of powder in plastic bottle
620628	100 g	100 g of powder in plastic bottle

### pH OF THE MEDIUM

7.2 ± 0.2

### USE

Chromatic™ *Bacillus cereus* is a chromogenic medium used with supplements for the isolation, differentiation and enumeration of *Bacillus cereus* from a wide variety of foods and environmental samples in the area of food production and handling

### TECHNIQUE

Refer to technical sheet of the product

### APPEARANCE OF THE MEDIUM

#### Powder medium

Appearance: free-flowing, homogeneous,

Colour: yellowish

#### Ready-to-use medium

Appearance: slightly opalescent

Colour: light straw

### SHELF LIFE

2 years

### QUALITY CONTROL

- Control of general characteristics, label and print
- Microbiological control  
Supplement: Chromatic™ *Bacillus cereus* Supplement  
Inoculum for productivity: 50-100 CFU  
Inoculum for selectivity: 10<sup>4</sup>-10<sup>6</sup> CFU  
Incubation Conditions: 18-24 h at 30 ± 2°C, in aerobiosis

Microorganism		Growth	Colony Appearance
<i>Bacillus cereus</i>	ATCC® 11778	Good	Blue-green colonies with opaque halo
<i>Bacillus subtilis</i>	ATCC® 6633	Inhibited	---
<i>Escherichia coli</i>	ATCC® 25922	Inhibited	---

### TABLE OF SYMBOLS

 Batch code	 Do not reuse	 Manufacturer	 Use by	 Fragile, handle with care
 Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult instructions for use	

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**LIOFILCHEM® S.r.l.**