CHROMATIC™ ESBL AGAR BASE

Chromogenic medium for detection of ESBLs and AmpC in Enterobacteriaceae directly from clinical specimen.

TYPICAL FORMULA	(g/l)
Peptone Mix	43.2
Chromogenic Mix	1.0
Agar	15.0
Final pH 7.2 ± 0.2 at 25°C	

DESCRIPTION

CHROMATIC™ ESBL AGAR BASE is a chromogenic medium used with supplements for detection of Extended-spectrum β -lactamase (ESBL) and AmpC producers.

ESBLs are enzymes that hydrolyze most penicillins and cephalosporins, inhibited by β -lactamase inhibitors such as clavulanic acid, sulbactam and tazobactam. AmpC-type enzymes hydrolyze penicillins, cephalosporins and monobactam and are poorly inhibited by the classical ESBL inhibitors, especially clavulanic acid.

PRINCIPLE

Peptones supply amino acids, nitrogen, carbon, minerals, vitamins and other nutrients which support the growth of microorganisms. Chromogenic mix allows the identification of microorganisms on the basis of the colony color and morphology. Agar is the solidifying agent.

The medium must be supplied with one of the following supplements:

- Chromatic™ ESBL Supplement (ref. 81089) contains a selective mix that inhibits the ESBL-non-producing organisms, including AmpC producers
- Chromatic™ ESBL+AmpC Supplement (ref. 81090) contains a selective mix that inhibits the ESBL-non-producing organisms, but allows the growth of AmpC producers.

PREPARATION

Suspend 59.2 g of powder in one liter of deionized or distilled water. Bring to boil and shake until completely dissolved. Sterilize at 121°C for 15 minutes. Cool up to 45-50°C. Aseptically, add rehydrated content of 2 vials (10 ml) of either Chromatic™ ESBL Supplement (ref. 81089) or Chromatic™ ESBL+AmpC Supplement (ref. 81090), as desired*. Pour in Petri dishes.

*Notice that Chromatic™ ESBL Supplement gives more selectivity compared to Chromatic™ ESBL+AmpC Supplement.

TECHNIQUE

Inoculate the plates by streaking directly the specimen onto the agar surface. Incubate aerobically at 37°C for 18-24 hours.

INTERPRETATION OF RESULTS

E. coli produces pink-reddish-mauve colonies.

Klebsiella spp, *Enterobacter* spp, *Serratia* spp produce green-blue colonies.

Proteus spp produces brown colonies.

STORAGE AND TRANSPORT CONDITIONS

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident. Store prepared plates at 2-8°C away from light.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *in vitro* diagnostic use only and must be used by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to the national and local regulations in force.

REFERENCES

- EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance. Version 1.0, 2013.
- Podschun R, Ullman U. *Klebsiella* spp as Nosocomial Pathogens: Epidemiology, Taxonomy, Typing Methods, and Pathogenicity Factors. *Clinical Microbiology Reviews*. 1998; 11 (4): 589–603.
- Geiss H.K. Comparison of two test kits for rapid identification of *Escherichia coli* by a beta-glucuronidase assay. *European Journal of Clinical Microbiology & Infectious Diseases*. 1990; 9 (2):151-152.

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10-30°C

PACKAGING

Ref.	Content	Packaging
610629	500 g	500 g of powder in plastic bottle
620629	100 g	100 g of powder in plastic bottle

pH OF THE MEDIUM

7.2 ± 0.2

USE

CHROMATIC™ ESBL AGAR BASE is a chromogenic medium used with supplements for detection of Extended-spectrum β-lactamase (ESBL) and AmpC producers

TECHNIQUE

Refer to technical sheet of the product

APPEARANCE OF THE MEDIUM

Powder medium

Appearance: fine, dry, homogeneous, free of extraneous material

Colour: beige

Ready-to-use medium

Appearance: slightly opalescent

Colour: amber

SHELF LIFE

2 years

QUALITY CONTROL

- Control of general characteristics, label and print
- Microbiological control
Supplement: Chromatic™ ESBL Supplement
Inoculum for productivity: 10-100 CFU/ml
Inoculum for selectivity: 10⁴-10⁵ CFU/ml
Incubation Conditions: 18-24 h at 35 ± 2°C, in aerobiosis

Microorganism

Microorganism		Growth	Colony colour
<i>Escherichia coli</i> (ESBL+, AmpC-)	DSM 22311	Good	Reddish
<i>Klebsiella pneumoniae</i> (ESBL+, AmpC-)	ATCC® 700603	Good	Green-blue
<i>Klebsiella pneumoniae</i> (ESBL-, AmpC+)	ATCC® BAA-1144	Inhibited	---
<i>Escherichia coli</i> (ESBL-, AmpC-)	ATCC® 25922	Inhibited	---

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult instructions for use	 Do not reuse

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1. EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance. Version 1.0, 2013.
2. Podschun R, Ullman U. *Klebsiella* spp as Nosocomial Pathogens: Epidemiology, Taxonomy, Typing Methods, and Pathogenicity Factors. *Clinical Microbiology Reviews*. 1998; 11 (4): 589-603.
3. Geiss H.K. Comparison of two test kits for rapid identification of *Escherichia coli* by a beta-glucuronidase assay. *European Journal of Clinical Microbiology & Infectious Diseases*. 1990; 9 (2):151-152.

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LIOFILCHEM® S.r.l.