

# Liofilchem® DIP-SLIDE

## Dip-Slide for the microbiological screening in **Urinary Tract Infections**



Description	Packaging	Ref.
<b>URITEST</b> (MAC CONKEY AGAR / CETRIMIDE AGAR / CLED ANDRADE AGAR) Enterobacteria, <i>Pseudomonas</i> spp., total bacterial count.	20 slides 120 slides	500152 51015
<b>URITEST N</b> (MAC CONKEY AGAR / CETRIMIDE AGAR / CLED AGAR) Enterobacteria, <i>Pseudomonas</i> spp., total bacterial count.	20 slides 120 slides 500 slides	500232 51023 51123
<b>URITEST C</b> (MAC CONKEY MUG AGAR / CETRIMIDE AGAR / CLED ANDRADE AGAR) Enterobacteria, <i>Escherichia coli</i> , <i>Pseudomonas</i> spp., total bacterial count.	20 slides 120 slides	500242 51024
<b>URITEST EC</b> (MAC CONKEY AGAR / E.COLI AGAR / CLED AGAR) Enterobacteria, <i>Escherichia coli</i> and coliforms, total bacterial count.	20 slides 120 slides	500412 51041
<b>URITEST 2</b> (CLED AGAR/ MAC CONKEY AGAR) Total bacterial count, enterobacteria.	20 slides 120 slides 500 slides	500302 51030 51130
<b>URITEST M</b> (MAC CONKEY AGAR / MALT SELECTIVE AGAR / CLED ANDRADE AGAR) Enterobacteria, yeasts, total bacterial count.	20 slides 120 slides 500 slides	500182 51018 51118
<b>URITEST MALTO</b> (MAC CONKEY AGAR / MALT AGAR / CLED AGAR) Enterobacteria, yeasts, total bacterial count.	20 slides 120 slides 500 slides	500402 51040 51140
<b>URITEST EF</b> (MAC CONKEY AGAR / BILE ESCULIN AGAR / CLED AGAR) Enterobacteria, enterococci, total bacterial count.	20 slides 120 slides	500702 51070
<b>URITEST PENTA</b> (CLED AGAR / MAC CONKEY AGAR / CETRIMIDE AGAR / MALT AGAR / BILE ESCULIN AGAR) Total bacterial count, enterobacteria, <i>Pseudomonas</i> spp., yeasts, enterococci.	20 slides 120 slides	500142 51014

## Dip-Slide for the microbiological screening in **Urogenital Infections**

Description	Packaging	Ref.
<b>VAGITEST</b> (CHOCOLATE / ROGOSA / THAYER MARTIN) Fastidious microorganisms, lactobacilli, <i>Neisseria gonorrhoeae</i> .	10 slides 120 slides	50020 51020

## Dip-Slide for the microbiological screening in **Skin Infections**

Description	Packaging	Ref.
<b>DERMATEST</b> (MALTO / CORN MEAL / DTM-TAPLIN) Pathogenic dermatophytes and others non pathogenic fungi.	10 slides 120 slides	50021 51021

Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.

[www.liofilchem.net](http://www.liofilchem.net)



**Liofilchem s.r.l.**  
Zona Industriale  
64026, Roseto degli Abruzzi (Te), Italy  
tel +39 0858930745 fax +39 0858930330

ref. 6550616  
Rev. 4 - 05.06.2014  
© Liofilchem 2014  
printed in Italy