

CultiControl™ freeze-dried microorganisms

Packaging: 1 vial containing 5 pellets

Non-enumerated CFU

Applications: Culture purposes, QC of ID devices, QC of AST devices



Description	CultiControl™ Ref.	notes
<i>Aggregatibacter aphrophilus</i> ATCC® 7901™*	89091	
<i>Aspergillus brasiliensis</i> ATCC® 16404™*	89021	
<i>Bacillus cereus</i> ATCC® 11778™*	89022	
<i>Bacillus subtilis</i> subsp. spizizenii ATCC® 6633™*	89023	
<i>Bacteroides fragilis</i> ATCC® 25285™*	89078	
<i>Bacteroides thetaiotaomicron</i> ATCC® 29741™*	89079	
<i>Campylobacter jejuni</i> subsp. jejuni ATCC® 33291™*	89086	
<i>Candida albicans</i> ATCC® 10231™*	89024	
<i>Candida albicans</i> ATCC® 90028™*	89072	
<i>Candida parapsilosis</i> ATCC® 22019™*	89071	
<i>Clostridium difficile</i> ATCC® 9689™*	89090	produces cytotoxin
<i>Clostridium perfringens</i> ATCC® 13124™*	89053	
<i>Clostridium sordellii</i> ATCC® 9714™*	89059	
<i>Enterobacter cloacae</i> subsp. cloacae ATCC® BAA-1143™*	89065	control strain for the AmpC disk test; strong positive
<i>Enterococcus faecalis</i> ATCC® 19433™*	89025	
<i>Enterococcus faecalis</i> ATCC® 29212™*	89026	
<i>Enterococcus faecalis</i> ATCC® 49532™*	89066	high level Gentamicin-resistant and Streptomycin-sensitive
<i>Enterococcus faecalis</i> ATCC® 49533™*	89067	high level Gentamicin-sensitive and Streptomycin-resistant
<i>Escherichia coli</i> ATCC® 25922™*	89027	
<i>Escherichia coli</i> ATCC® 8739™*	89028	
<i>Escherichia coli</i> NCTC 11954	89068	beta lactamase producer
<i>Haemophilus influenzae</i> ATCC® 49766™*	89076	
<i>Haemophilus influenzae</i> ATCC® 49247™*	89077	
<i>Issatchenkia orientalis</i> ATCC® 6258™*	89073	
<i>Klebsiella pneumoniae</i> ATCC® BAA-1705™*	89088	Modified Hodge Test (MHT) positive control
<i>Klebsiella pneumoniae</i> ATCC® BAA-1706™*	89087	Modified Hodge Test (MHT) negative control
<i>Klebsiella pneumoniae</i> ATCC® BAA-2146™*	89069	New Delhi metallo-beta-lactamase (NDM-1) positive
<i>Klebsiella pneumoniae</i> subsp. pneumoniae ATCC® 700603™*	89070	ESBL positive
<i>Lactobacillus acidophilus</i> ATCC® 4356™*	89080	
<i>Lactobacillus paracasei</i> subsp. paracasei ATCC® BAA-52™*	89055	
<i>Lactococcus lactis</i> subsp. lactis ATCC® 19435™*	89082	
<i>Lactobacillus leichmannii</i> ATCC® 4797™*	89081	
<i>Legionella pneumophila</i> subsp. pneumophila ATCC® 33152™*	89052	
<i>Klebsiella pneumoniae</i> subsp. pneumoniae ATCC® 13883™*	89089	
<i>Listeria innocua</i> ATCC® 33090™*	89029	
<i>Listeria ivanovii</i> subsp. ivanovii ATCC® 19119™*	89030	
<i>Listeria monocytogenes</i> ATCC® 19111™*	89031	serotype 1
<i>Listeria monocytogenes</i> ATCC® 19115™*	89051	serotype 4b
<i>Listeria monocytogenes</i> ATCC® 13932™*	89085	serotype 4b
<i>Listeria monocytogenes</i> ATCC® 7644™*	89060	
<i>Neisseria gonorrhoeae</i> ATCC® 19424™*	89074	



* The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Liofilchem® s.r.l. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

The ATCC Licensed Derivative® program was created to ensure quality and address global responsibilities for handling biological materials. A royalty fee will be charged on behalf of ATCC® for products derived from ATCC® biological material. In order to receive these products, an authorized representative from your organization must agree to the terms of the End User Agreement. It is mandatory that your organization registers for this program before we are able to process your order for lyophilized microorganism preparations containing ATCC® biological material. If you have additional questions about the program contact Liofilchem Customer Service.

CultiControl™ freeze-dried microorganisms

Packaging: 1 vial containing 5 pellets

Non-enumerated CFU

Applications: Culture purposes, QC of ID devices, QC of AST devices



Description	CultiControl™ Ref.	notes
<i>Neisseria gonorrhoeae</i> ATCC® 31426™*	89075	beta lactamase producer
<i>Plesiomonas shigelloides</i> ATCC® 14029™*	89094	
<i>Proteus mirabilis</i> ATCC® 25933™*	89032	
<i>Proteus mirabilis</i> ATCC® 12453™*	89049	
<i>Proteus mirabilis</i> ATCC® 29906™*	89083	
<i>Pseudomonas aeruginosa</i> ATCC® 27853™*	89033	
<i>Pseudomonas aeruginosa</i> ATCC® 9027™*	89034	
<i>Rhodococcus equi</i> ATCC® 6939™*	89035	recommended for CAMP test for <i>Listeria monocytogenes</i>
<i>Saccharomyces cerevisiae</i> ATCC® 9763™*	89036	
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Enteritidis ATCC® 13076™*	89084	group D
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium ATCC® 14028™*	89037	
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium ATCC® 13311™*	89054	
<i>Shigella flexneri</i> ATCC® 12022™*	89038	serotype 2b
<i>Staphylococcus aureus</i> ATCC® 33862™*	89042	recommended for CAMP test
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> ATCC® 29213™*	89041	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> ATCC® 25923™*	89040	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> ATCC® 43300™*	89043	methicillin resistant
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> ATCC® 700699™*	89093	Methicillin resistant; Mu50; reduced Vancomycin susceptibility
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> ATCC® 700698™*	89092	Methicillin resistant; GRD MIC Test Strip control
<i>Staphylococcus aureus</i> ATCC® 6538™*	89044	
<i>Staphylococcus aureus</i> NCTC 12493	89039	
<i>Staphylococcus epidermidis</i> ATCC® 12228™*	89045	
<i>Streptococcus agalactiae</i> ATCC® 13813™*	89046	group B; non-hemolytic in absence of CAMP Factor
<i>Streptococcus bovis</i> ATCC® 33317™*	89061	
<i>Streptococcus mutans</i> ATCC® 25175™*	89062	
<i>Streptococcus pneumoniae</i> ATCC® 27336™*	89063	
<i>Streptococcus pneumoniae</i> ATCC® 49619™*	89047	low level penicillin resistance by oxacillin test
<i>Streptococcus pyogenes</i> ATCC® 19615™*	89048	group A
<i>Streptococcus sanguinis</i> ATCC® 10556™*	89064	
<i>Vibrio parahaemolyticus</i> ATCC® 17802™*	89056	
<i>Yersinia enterocolitica</i> subsp. <i>enterocolitica</i> ATCC® 9610™*	89050	biovar 1; serogroup O:8



* The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Liofilchem® s.r.l. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

The ATCC Licensed Derivative® program was created to ensure quality and address global responsibilities for handling biological materials. A royalty fee will be charged on behalf of ATCC® for products derived from ATCC® biological material. In order to receive these products, an authorized representative from your organization must agree to the terms of the End User Agreement. It is mandatory that your organization registers for this program before we are able to process your order for lyophilized microorganism preparations containing ATCC® biological material. If you have additional questions about the program contact Liofilchem Customer Service.

Liofilchem® CultiControl™

Quanti-CultiControl™ freeze-dried microorganisms

Packaging: 1 vial containing 1 pellet + 1 rehydration fluid vial

Quantitative: <100 CFU / 0.1 mL inoculum

Applications: Growth Promotion Testing, Microbial Limits Testing, Microbial Enumeration Testing

Description	Quanti-CultiControl™ Ref.	notes
<i>Aspergillus brasiliensis</i> ATCC® 16404™*	89501	
<i>Bacillus cereus</i> ATCC® 10876™*	89502	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> ATCC® 6633™*	89503	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NCIMB 8054	89504	
<i>Bacteroides fragilis</i> ATCC® 25285™*	89505	
<i>Brevundimonas diminuta</i> ATCC® 19146™*	89506	
<i>Burkholderia cepacia</i> ATCC® 25416™*	89507	
<i>Candida albicans</i> ATCC® 10231™*	89508	
<i>Candida albicans</i> ATCC® 14053™*	89509	
<i>Candida albicans</i> ATCC® 2091™*	89510	
<i>Candida albicans</i> NCYC 1363	89511	
<i>Clostridium perfringens</i> ATCC® 13124™*	89512	
<i>Clostridium sporogenes</i> ATCC® 11437™*	89513	
<i>Clostridium sporogenes</i> ATCC® 19404™*	89514	
<i>Clostridium sporogenes</i> NCIMB 12343	89515	
<i>Enterobacter aerogenes</i> ATCC® 13048™*	89516	
<i>Enterococcus faecalis</i> ATCC® 29212™*	89517	
<i>Enterococcus faecalis</i> ATCC® 7080™*	89518	
<i>Escherichia coli</i> ATCC® 8739™*	89519	
<i>Geobacillus stearothermophilus</i> ATCC® 12980™*	89521	
<i>Geobacillus stearothermophilus</i> ATCC® 7953™*	89522	
<i>Kocuria rhizophila</i> ATCC® 9341™*	89523	
<i>Lactobacillus fermentum</i> ATCC® 9338™*	89524	
<i>Listeria monocytogenes</i> ATCC® 19115™*	89525	serotype 4b
<i>Micrococcus luteus</i> ATCC® 4698™*	89526	
<i>Pseudomonas aeruginosa</i> ATCC® 27853™*	89527	
<i>Pseudomonas aeruginosa</i> ATCC® 9027™*	89528	
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Choleraesuis ATCC® 10708™*	89529	
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium ATCC® 14028™*	89531	
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium ATCC® 13311™*	89530	
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Abony NCTC 6017	89532	serovar Abony
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> ATCC® 25923™*	89533	
<i>Staphylococcus aureus</i> ATCC® 6538™*	89535	
<i>Staphylococcus aureus</i> ATCC® 6538P™*	89534	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> NCIMB 9518	89536	
<i>Staphylococcus epidermidis</i> ATCC® 12228™*	89537	
<i>Streptococcus pyogenes</i> ATCC® 19615™*	89538	group A



* The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Liofilchem® s.r.l. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

The ATCC Licensed Derivative® program was created to ensure quality and address global responsibilities for handling biological materials. A royalty fee will be charged on behalf of ATCC® for products derived from ATCC® biological material. In order to receive these products, an authorized representative from your organization must agree to the terms of the End User Agreement. It is mandatory that your organization registers for this program before we are able to process your order for lyophilized microorganism preparations containing ATCC® biological material. If you have additional questions about the program contact Liofilchem Customer Service.