


Liofilchem®

CONTACT SLIDE

flexible contact slide for microbiological screening of surfaces

ISO 18593



Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.

Description	Packaging	Ref.	# Product with minimum order
CONTACT SLIDE Chrom 1 (Chromatic™ E.coli-coliform / Chromatic™ Salmonella) Chromogenic selective media for detecting <i>E. coli</i> , coliforms and <i>Salmonella</i> spp. on surfaces or in food and waters.	20 slides 120 slides	525292 # 53529 #	
CONTACT SLIDE Chrom 2 (Chromatic™ E.coli-coliform / PCA+TTC + NEUTRALIZING) Chromogenic selective medium for detecting <i>E. coli</i> , coliforms and for total bacterial count.	20 slides 120 slides	525392 53539	
CONTACT SLIDE Chrom 3 (Chromatic™ E.coli-coliform / BAIRD PARKER AGAR + NEUTRALIZING) Chromogenic selective medium for detecting <i>E. coli</i> , coliforms and <i>Staphylococcus aureus</i> with inactivation of disinfectants.	20 slides 120 slides	525462 # 53546 #	
CONTACT SLIDE 1 (PCA + TTC + NEUTRALIZING / V.R.B.G. AGAR + NEUTRALIZING) Total bacterial count, enterobacteria, with inactivation of disinfectants.	20 slides 120 slides	525262 53526	
CONTACT SLIDE 2 (PCA + TTC + NEUTRALIZING / ROSE BENGAL CAF AGAR + NEUTRALIZING) Total bacterial count, yeasts and moulds, with inactivation of disinfectants.	20 slides 120 slides	525272 53527	
CONTACT SLIDE 3 (PCA+TTC + NEUTRALIZING / VOGEL JOHNSON AGAR) Total bacterial count with inactivation of disinfectants, staphylococci.	20 slides 120 slides	525302 53530	
CONTACT SLIDE 4 (CETRIMIDE AGAR + NEUTRALIZING / ROSE BENGAL CAF AGAR + NEUTRALIZING) <i>Pseudomonas</i> spp., yeasts and moulds.	20 slides 120 slides	525282 # 53528 #	
CONTACT SLIDE 5 (V.R.B.G. AGAR + NEUTRALIZING / BILE ESCULIN AGAR + NEUTRALIZING) Enterobacteria, fecal streptococci, with inactivation of disinfectants.	20 slides 120 slides	525312 53531	
CONTACT SLIDE 8 (PCA + TTC + NEUTRALIZING / V.R.B.L. AGAR + NEUTRALIZING) Total bacterial count, coliforms, with inactivation of disinfectants.	20 slides 120 slides	525622 53562	
CONTACT SLIDE 9 (PCA + TTC + NEUTRALIZING / PCA + TTC + NEUTRALIZING) Total bacterial count with inactivation of disinfectants.	20 slides 120 slides	525322 53532	
CONTACT SLIDE 10 (ROSE BENGAL CAF AGAR + NEUTRALIZING / VOGEL JOHNSON AGAR) Yeasts and moulds with inactivation of disinfectants, staphylococci.	20 slides 120 slides	525452 # 53545 #	
CONTACT SLIDE TSA + NEUTRALIZING / TSA + NEUTRALIZING Total bacterial count with inactivation of disinfectants.	20 slides 120 slides	525472 # 53547 #	
CONTACT SLIDE TSA + NEUTRALIZING / ROSE BENGAL CAF AGAR + NEUTRALIZING Total bacterial count with inactivation of disinfectants, yeasts and moulds.	20 slides 120 slides	525482 # 53548 #	



Liofilchem s.r.l.
Zona Industriale
64026, Roseto degli Abruzzi (Te), Italy
tel +39 0858930745 fax +39 0858930330

www.liofilchem.net

ref. 6550616
Rev. 5.1 - 30.06.2014
© Liofilchem 2014
printed in Italy

Liofilchem® DIP-SLIDE

Dip-Slide for the microbiological screening in **Urinary Tract Infections**



Description	Packaging	Ref.
URITEST (MAC CONKEY AGAR / CETRIMIDE AGAR / CLED ANDRADE AGAR) Enterobacteria, <i>Pseudomonas</i> spp., total bacterial count.	20 slides 120 slides	500152 51015
URITEST N (MAC CONKEY AGAR / CETRIMIDE AGAR / CLED AGAR) Enterobacteria, <i>Pseudomonas</i> spp., total bacterial count.	20 slides 120 slides 500 slides	500232 51023 51123
URITEST C (MAC CONKEY MUG AGAR / CETRIMIDE AGAR / CLED ANDRADE AGAR) Enterobacteria, <i>Escherichia coli</i> , <i>Pseudomonas</i> spp., total bacterial count.	20 slides 120 slides	500242 51024
URITEST EC (MAC CONKEY AGAR / E.COLI AGAR / CLED AGAR) Enterobacteria, <i>Escherichia coli</i> and coliforms, total bacterial count.	20 slides 120 slides	500412 51041
URITEST 2 (CLED AGAR/ MAC CONKEY AGAR) Total bacterial count, enterobacteria.	20 slides 120 slides 500 slides	500302 51030 51130
URITEST M (MAC CONKEY AGAR / MALT SELECTIVE AGAR / CLED ANDRADE AGAR) Enterobacteria, yeasts, total bacterial count.	20 slides 120 slides 500 slides	500182 51018 51118
URITEST MALTO (MAC CONKEY AGAR / MALT AGAR / CLED AGAR) Enterobacteria, yeasts, total bacterial count.	20 slides 120 slides 500 slides	500402 51040 51140
URITEST EF (MAC CONKEY AGAR / BILE ESCULIN AGAR / CLED AGAR) Enterobacteria, enterococci, total bacterial count.	20 slides 120 slides	500702 51070
URITEST PENTA (CLED AGAR / MAC CONKEY AGAR / CETRIMIDE AGAR / MALT AGAR / BILE ESCULIN AGAR) Total bacterial count, enterobacteria, <i>Pseudomonas</i> spp., yeasts, enterococci.	20 slides 120 slides	500142 51014

Dip-Slide for the microbiological screening in **Urogenital Infections**

Description	Packaging	Ref.
VAGITEST (CHOCOLATE / ROGOSA / THAYER MARTIN) Fastidious microorganisms, lactobacilli, <i>Neisseria gonorrhoeae</i> .	10 slides 120 slides	50020 51020

Dip-Slide for the microbiological screening in **Skin Infections**

Description	Packaging	Ref.
DERMATEST (MALTO / CORN MEAL / DTM-TAPLIN) Pathogenic dermatophytes and others non pathogenic fungi.	10 slides 120 slides	50021 51021

Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.

www.liofilchem.net



Liofilchem s.r.l.
Zona Industriale
64026, Roseto degli Abruzzi (Te), Italy
tel +39 0858930745 fax +39 0858930330

ref. 6550616
Rev. 4 - 05.06.2014
© Liofilchem 2014
printed in Italy